

## QC Wireless SDK Instructions for Use

update note:

### 1.introduce

#### 1.1The role of the SDK

### 2. API description

#### 2.0 Access Conditions

#### 2.1 Permissions required by the SDK

#### 2.2 Access conditions

#### 2.3.1API

Scanning Devices

device connection:BleOperateManager.getInstance()

#### 2.3.2 Feature list:

Synchronize time, get the list of functions supported by the device

bracelet battery

Continuous heart rate, blood oxygen, blood pressure switch

Set watch sports goals

Find equipment

R11 RING Display settings

Ring photo control

Set the Ring to factory reset

Shutdown the Ring

Wearing Calibration

Close Bluetooth

#### 2.3.3 Data synchronization:

Synchronize steps, distance, kcal for the day

Synchronized step data details

Synchronize the details of new sleep data and return according to SetTimeRsp

Sedentary data synchronization

Sync heart rate data

Synchronized blood pressure function

Synchronized blood oxygen function

Synchronized Time Interval blood oxygen function

Synchronized Time Interval temperature function

Synchronized pressure function

Synchronized hrv function

Synchronized skin(body) temperature function

Synchronous training records

#### 2.3.7Manual measurement

#### 2.3.8 Touch and gestures

#### 2.3.9 Changes in Ring measurement data are proactively reported to the APP

#### 2.3.10 APP opens exercise type

Muslim Data Synchronization

Setting Ring user Id

设置个人信息

## QC Wireless SDK Instructions for Use

---

1. Author: James

update note:

- 1. (2021/07/06) scan, connect, measure commands
- 2. (2021/07/20) Add setting command
- 3. (2021/07/21) Increase step count, heart rate, sleep data sync
- 4. (2022/02/28) Add new sleep algorithm
- 5. (2023/03/10) Add message switch synchronization, body temperature, user information, watch manual measurement of heart rate, blood pressure results
- 6. (2023/03/16) Add bt connection and contact person
- 7. (2024/02/23) Add manual pressure, APP movement, pressure synchronization, pressure setting
- 8. (2025/05/11) Add One-click measurement
- 9. (2025/05/11) Temperature switch setting

1.introduce

1.1The role of the SDK

Provide partner companies with the Android Bluetooth SDK for use with Green Orange wireless devices that provide basic and advanced functionality for a major watch or other device. This document is intended to explain the usage context, functionality, etc. of the API. Intended Audience and Reading Recommendations The intended audience and reader recommendations in this article are shown in Annex 1.

Reader	Role
Software Architecture Engineer	Architecture Analysis and Technical Guidance
Android development engineer	Have a certain android development ability, understand Ble related development technology

2. API description

2.0 Access Conditions

Android 5.0 or above, Bluetooth 4.0 or above.

2.1 Permissions required by the SDK

```
//network permissions
<uses-permission android:name="android.permission.INTERNET" />
//Bluetooth related permissions
<uses-permission android:name="android.permission.BLUETOOTH" />
<uses-permission android:name="android.permission.BLUETOOTH_ADMIN" />
//Storage related permissions
<uses-permission android:name="android.permission.WRITE_EXTERNAL_STORAGE" />
<uses-permission android:name="android.permission.READ_EXTERNAL_STORAGE" />
//Location permission
<uses-permission android:name="android.permission.ACCESS_COARSE_LOCATION" />
<uses-permission android:name="android.permission.ACCESS_FINE_LOCATION" />
```

If it is android 12 or higher system

```
<uses-permission android:name="android.permission.BLUETOOTH_CONNECT" />
<uses-permission android:name="android.permission.BLUETOOTH_SCAN" />
<uses-permission android:name="android.permission.BLUETOOTH_ADVERTISE" />
```

## 2.2 Access conditions

- Green Orange Wireless Wearables
- Green Orange Wireless SDK and documentation

### 2.3.1 API

#### Scanning Devices

```
//start scan
BleScannerHelper.getInstance().scanDevice(final Context context, UUID mUuid,
final ScanWrapperCallback scanCallback);
//stop scan
BleScannerHelper.getInstance().stopScan(Context context)
//Specified device scan
BleScannerHelper.getInstance().scanTheDevice(final Context context, final
String macAddress, final OnTheScanResult scanResult)
```

#### device connection:BleOperateManager.getInstance()

```
//direct connection
BleOperateManager.getInstance().connectDirectly(smartwatch.deviceAddress)
//scan connections
BleOperateManager.getInstance().connectWithScan(smartwatch.deviceAddress)
//disconnect
BleOperateManager.getInstance().unBindDevice()
//reconnect
BleOperateManager.getInstance().setNeedConnect(boolean needConnect)
//Called when bluetooth is turned off
BleOperateManager.getInstance().setBluetoothTurnOff(false)
BleOperateManager.getInstance().disconnect()
//Turn on the system bluetooth monitor
BleOperateManager.getInstance().setBluetoothTurnOff(true)

//bind seccess
CommandHandle.getInstance()
    .executeReqCmdNoCallback(
        SimpleKeyReq(Constants.CMD_BIND_SUCCESS)
    )
```

## 2.3.2 Feature list:

### Synchronize time, get the list of functions supported by the device

```
//set time
CommandHandle.getInstance().executeReqCmd(SetTimeReq(0),
ICommandResponse<SetTimeRsp>() {}))
//Callback Description
public class SetTimeRsp extends BaseRspCmd {
    //support body temperature
    public boolean mSupportTemperature;
    //watch face
    public boolean mSupportPlate;
    //Support the menstrual cycle
    public boolean mSupportMenstruation;
    //Support custom watch faces
    public boolean mSupportCustomWallpaper;
    //Support blood oxygen
    public boolean mSupportBloodOxygen;
    //blood pressure support
    public boolean mSupportBloodPressure;
    //Support fatigue
    public boolean mSupportFeature;
    //Support one-key detection
    public boolean mSupportOneKeyCheck;
    //weather support
    public boolean mSupportWeather;
    //Support for new sleep protocol
    public boolean mNewSleepProtocol;
    //Supports up to 6 or 3 dials
    public int mMaxWatchFace;
    // hrv support
    public boolean mSupportHrv;
}
```

Functional support

```
CommandHandle.getInstance()
    .executeReqCmd(
        DeviceSupportReq.getReadInstance(),
        ICommandResponse<DeviceSupportFunctionRsp> {

    //support touch
    public boolean supportTouch;
    //support muslim
    public boolean supportMoslin;
    //support ble pair
    public boolean supportAPPRvision;
    //support Heart rate calibration
    public boolean supportBlePair;
    //support gesture
    public boolean supportGesture;
    //support music
    public boolean supportRingMusic;
    //support video
    public boolean supportRingVideo;
```

```

//support ebook
public boolean supportRingEbook;
//support camera
public boolean supportRingCamera;
//support phone call
public boolean supportRingPhoneCall;
//support game
public boolean supportRingGame;
}

```

## bracelet battery

```

        CommandHandle.getInstance().executeReqCmd(new
SimpleKeyReq(Constants.CMD_GET_DEVICE_ELECTRICITY_VALUE), new
ICommandResponse<BatteryRsp>() {

@Override
public void onDataResponse(BatteryRsp resultEntity) {
    if (resultEntity.getStatus() == BaseRspCmd.RESULT_OK) {
        //Get battery successfully
    }
}

});

//Callback Description
public class BatteryRsp extends BaseRspCmd {
    //battery 【0-100】
    private int batteryValue;
}

```

## Continuous heart rate, blood oxygen, blood pressure switch

```

//Read continuous heart rate settings (too Low/High Reminder is notify in
DeviceNotifyListener 0x3A.need
BleOperateManager.getInstance().addOutDeviceListener(ListenerKey.Heart,yourListen
er))

CommandHandle.getInstance()
    .executeReqCmd(HeartRateSettingReq.getReadInstance(),
        ICommandResponse<HeartRateSettingRsp> {
//it.isEnabled        switch
//it.heartInterval    heart rate measurement interval
//it.startInterval    minValue/minUnit 5、10、30
//tooLowReminder low value reminder 0、40、45、50    0:tooLowReminder is off
other value: is on
//tooHighReminder high value reminder 0、110、120、130、140、150
0:tooHighReminder is off    other value: is on
// mainSwitch 0:no support main control 1:support main control and is on
2:1:support main control and is off
})

//Read Continuous SpO2 settings
CommandHandle.getInstance()
    .executeReqCmd(BloodOxygenSettingReq.getReadInstance(),

```

```

        ICommandResponse<BloodOxygenSettingRsp> {
//it.isEnabled      switch
//it.interval      interval time unit min
    })

    //Read continuous blood pressure settings
    CommandHandle.getInstance()
        .executeReqCmd(BpSettingReq.getReadInstance(),
            ICommandResponse<BpSettingRsp> {

    })

    //Read pressure setting
    CommandHandle.getInstance()
        .executeReqCmd(PressureSettingReq.getReadInstance(),
            ICommandResponse<PressureSettingRsp>() {
switch:it.isEnabled
        })

        //Write continuous heart rate switch isEnabled: true on,
false: off  hrInterval: 10, 15, 20, 30, 60

        CommandHandle.getInstance().executeReqCmd(
            HeartRateSettingReq.getWriteInstance(true,hrInterval),
            ICommandResponse<HeartRateSettingRsp> {

    })

    //Write continuous heart rate switch isEnabled: true on, false: off
hrInterval: 10, 15, 20, 30, 60
    //startInterval minValue/minUnit 5、10、30
    //tooLowReminder low value reminder 0、40、45、50 0:tooLowReminder is
off  other value: is on
    //tooHighReminder high value reminder 0、110、120、130、140、150
0:tooHighReminder is off  other value: is on
    // mainSwitch 0:no support main control 1:support main control and is on
2:support main control and is off
    // too Low/High Reminder is notify in DeviceNotifyListener 0x3A(need
BLEOperateManager.getInstance().addOutDeviceListener(ListenerKey.Heart,yourListen
er)

    CommandHandle.getInstance().executeReqCmd(
        HeartRateSettingReq.getWriteInstance(
            isEnabled,
            hrInterval,
            startInterval,
            tooLowReminder,
            tooHighReminder
        ),

            ICommandResponse<HeartRateSettingRsp> {

    })

    //add Heart rate tooLow/High Reminder Listener
    BLEOperateManager.getInstance()
        .addOutDeviceListener(ListenerKey.Heart,
myDeviceNotifyListener)

    //Write continuous blood oxygen switch isEnabled: true on, false: off
    CommandHandle.getInstance().executeReqCmd(
        BloodOxygenSettingReq.getWriteInstance(boolean isEnabled),

```

```

        ICommandResponse<BloodOxygenSettingRsp> {

    })

    //write continuous interval blood oxygen switch isEnabled: true on, false:
off. interval: interval time unit min
    CommandHandle.getInstance().executeReqCmd(
        BloodOxygenSettingReq.getWriteInstance(boolean isEnabled,int interval),
        ICommandResponse<BloodOxygenSettingRsp> {

    })

    //write blood pressure switch

    CommandHandle.getInstance().executeReqCmd(BpSettingReq.getWriteInstance(boolean
isEnabled, StartEndTimeEntity startEndTimeEntity, int multiple),
        ICommandResponse<BpSettingRsp> {

    })
    BpSettingRsp,Parameter Description
isEnabled: true on false off
StartEndTimeEntity The parameter description is the same as above
multiple default 60

    //Write pressure setting switch,switch isEnabled: true on, false: off
    CommandHandle.getInstance().executeReqCmd(
        PressureSettingReq.getWriteInstance(isEnabled),
        ICommandResponse<PressureSettingRsp> {

    })

    //hrv readme
    CommandHandle.getInstance()
        .executeReqCmd(HrvSettingReq.getReadInstance(),
            ICommandResponse<HRVSettingRsp>() {

    })

    //hrv write
    CommandHandle.getInstance().executeReqCmd(
        HrvSettingReq(true),
        ICommandResponse<HRVSettingRsp> {

    })

    //hrv write interval
    CommandHandle.getInstance().executeReqCmd(
        HrvSettingReq(true,int interval),
        ICommandResponse<HRVSettingRsp> {

    })

    //body temperature write
    CommandHandle.getInstance().executeReqCmd(
        SugarLipidsSettingReq.getWriteInstance(0x03,isEnable,interval),
        ICommandResponse<BloodSugarLipidsSettingRsp> { })
    // temperature interval 2-120s
    CommandHandle.getInstance().executeReqCmd(
        SugarLipidsSettingReq.getWriteInstance(0x03,isEnable,0,interval),

```

```

        ICommandResponse<BloodSugarLipidsSettingRsp> { })

body temperature read
    CommandHandle.getInstance()
        .executeReqCmd(SugarLipidsSettingReq.getReadInstance(0x03),
            ICommandResponse<BloodSugarLipidsSettingRsp>() {
                if (it.type == 0x03.toByte()) {
                    //it.isEnabled
                }
            })

        //ring bluetooth enable 0:disable 1:enable
        CommandHandle.getInstance()
            .executeReqCmd(BluetoothEnableReq.getWriteInstance(int enable)

```

## Set watch sports goals

```

        CommandHandle.getInstance().executeReqCmd(
            TargetSettingReq.getWriteInstance(
                final int step, final int calorie, final int distance, final int
sportMinute, final int sleepMinute
            ), null
        )

```

### Parameter Description

step: step target

calorie: Calorie target target, unit kcal, write kcal to a\*1000

distance: Distance to target, in meters

sportMinute: Exercise duration target, in minutes

sleepMinute: Sleep duration target, in minutes

## Find equipment

```

CommandHandle.getInstance().executeReqCmd(FindDeviceReq(), null)

```

## R11 RING Display settings

```

// isEnabled: whether to display the switch
// isLeft: default left-hand mode (true = left hand, false = right hand)
// light: brightness level, default is 1
// maxLight: maximum brightness, default is 5
// start / end: screen-on and screen-off time
//             Full day: start = 0, end = 0
//             Otherwise, time is the number of minutes since midnight
//             For example, 3:00 AM = 180

CommandHandle.getInstance().executeReqCmd(palmScreen, ICommandResponse<PalmScreenR
sp>{

    })

//read

```



```
CommandHandle.getInstance()  
    .executeReqCmd(  
        PalmScreenReq.getRingReadInstance(),  
        ICommandResponse<PalmScreenRsp>() {  
  
            })
```

PalmScreenRsp description:

```
/**  
 * Whether the feature is enabled  
 */  
private boolean isEnabled;  
  
/**  
 * Whether the device is worn on the left hand  
 */  
private boolean isLeft;  
  
/**  
 * Whether palm gesture (wrist flip) is required  
 */  
private boolean needPalm;  
  
/**  
 * Screen brightness level  
 */  
private int screenLight;  
  
/**  
 * Maximum brightness level  
 */  
private int maxLight;  
  
/**  
 * All-day mode (Do Not Disturb)  
 */  
private boolean dnd;  
  
/**  
 * Start time  
 */  
private int start;  
  
/**  
 * End time  
 */  
private int end;
```

## Ring photo control

```
//The bracelet enters the camera interface

CommandHandle.getInstance().executeReqCmd(CameraReq(CameraReq.ACTION_INT0_CAMARA_
UI), null)

//The wristband is controlled by the bright screen on the camera
interface. The APP will send the bright screen command to keep the watch bright.
It is recommended to send it every 2 seconds.
CommandHandle.getInstance().executeReqCmd(
    CameraReq(CameraReq.ACTION_KEEP_SCREEN_ON),
    null
)
//Bracelet click to take a photo event
monitoring

BleOperateManager.getInstance().addNotifyListener(Constants.CMD_TAKING_PICTURE, new ICommandResponse<CameraNotifyRsp>(){

    @Override
    public void onDataResponse(CameraNotifyRsp resultEntity) {

    }

});

resultEntity.getAction()
Parameter Description
//The watch exits the camera interface
CameraNotifyRsp.ACTION_FINISH
//The watch clicked on the photo event
CameraNotifyRsp.ACTION_TAKE_PHOTO

CommandHandle.getInstance().executeReqCmd(CameraReq(CameraReq.ACTION_FINISH),
null)
```

## Set the Ring to factory reset

```
CommandHandle.getInstance().executeReqCmd(RestoreKeyReq(Constants.CMD_RE_STORE), null)
```

- message push

```
//The watch message push switch is fully turned on, and the APP should be
actively opened
CommandHandle.getInstance().executeReqCmd(
    SetANCSReq(), null
)
//Send message push to watch
MessPushUtil.pushMsg(type,message:String)

PushMsgUintReq parameter description
```

type:

0x00: Call reminder 0x01: SMS reminder 0x02: QQ reminder 0x03: WeChat reminder,  
0x04: incoming call to answer or hang up 0x05: Facebook message reminder  
0x06: WhatsApp message reminder  
0x07: Twitter message reminder 0x08: Skype message reminder 0x09: Line message reminder 0x0a: LinkedIn  
0x0b: Instagram 0x0c: TIM message 0x0d: Snapchat  
0x0e: others other types of notifications

## Shutdown the Ring

```
BleOperateManager.getInstance().shutdown()
```

## Wearing Calibration

```
start :enable:true
end :enable:false
    BleOperateManager.getInstance().ringCalibration(false
) {
    //2:Measuring 1:success 3:fail
    it.success
}

...

#### firmware version number, hardware version number
```java
//hardware information

CommandHandle.getInstance().execReadCmd(CommandHandle.getInstance().getReadHwRequest());

//firmware information

CommandHandle.getInstance().execReadCmd(CommandHandle.getInstance().getReadFmRequest());

Receiving implements this QCBluetoothCallbackCloneReceiver refer to demo
MyBluetoothReceiver
Judging UUID in the callback onCharacteristicRead

override fun onCharacteristicRead(uuid: String?, data: ByteArray?) {
    if (uuid != null && data != null) {
        val version = String(data, Charsets.UTF_8)
        when(uuid){
            Constants.CHAR_FIRMWARE_REVISION.toString() -> {
                //Firmware version number version
            }
            Constants.CHAR_HW_REVISION.toString() -> {
                //hardware version number version number version
            }
        }
    }
}
```

```
}
```

## Close Bluetooth

```
//Generally, charging will restart the Bluetooth.
```java
    CommandHandle.getInstance()
        .executeReqCmd(
            BluetoothCloseReq.getWriteInstance(),
object : ICommandResponse<BluetoothCloseRsp> {
    override fun onDataResponse(blueetoothCloseRsp: BluetoothCloseRsp) {
        //nothing
    }
})
```

### 2.3.3 Data synchronization:

#### Synchronize steps, distance, kcal for the day

```
CommandHandle.getInstance().executeReqCmd(
    SimpleKeyReq(Constants.CMD_GET_STEP_TODAY),
    ICommandResponse<TodaySportDataRsp> {})

TodaySportDataRsp parameter description
// days ago
private int daysAgo;
// date: year
private int year;
// date: month
private int month;
// date: day
private int day;
// total steps
private int totalSteps;
// running steps/aerobic steps
private int runningSteps;
// calorie value
private int calorie;
// walking distance
private int walkDistance;
// Movement time, in seconds
private int sportDuration;
// sleep time in seconds
private int sleepDuration;
```

#### Synchronized step data details

```
dayOffset 0: Today 1: Yesterday 2: The day before yesterday, supports
synchronization for up to 7 days
CommandHandle.getInstance().executeReqCmd(
    ReadDetailSportDataReq(dayOffset, 0, 95),
```

```

        ICommandResponse<ReadDetailSportDataRsp> {

    })
    BleStepDetails parameter description
    //year
    private int year;
    //moon
    private int month;
    //day
    private int day;
    //15 minutes a point, the total number of points in a day is 96 points, [0, 95],
    used to calculate the details of the number of steps per hour
    private int timeIndex=0;
    // calorie unit card
    private int calorie=0;
    //Step count
    private int walkSteps=0;
    //distance in meters
    private int distance=0;
    // keep for now
    private int runSteps=0;

```

## Synchronize the details of new sleep data and return according to SetTimeRsp

```

    //offset 0 today 1 yesterday
    public void syncSleepList(int offset, final ILargeDataSleepResponse response)
    SleepNewProtoResp parameter description
    //sleep start time
    private int st;
    //sleep end time
    private int et;
    //sleep collection
    private List<DetailBean> list;
    DetailBean parameter description
    // duration of a sleep type
    private int d;
    //0未睡眠, 1摘下, 2浅睡, 3深睡, 4眼动, 5清醒;
    private int t;

    //support lunch
    public void syncSleepList(int offset, final ILargeDataSleepResponse
    response, final ILargeDataLaunchSleepResponse lunchSleepResponse){

    }
    SleepNewProtoResp
    lunch start time
    private int lunchSt;
    lunch end time
    private int lunchEt;
    private boolean lunchBreak;
    private List<DetailBean> list;

```

## Sedentary data synchronization

```
//offset 0 today 1 yesterday
public void syncLongSitList(int offset, final ILargeSettingForLongDataResponse
response)
LongSitResp    parameter description
// offset 0 today 1 yesterday
private int index;
//Sedentary collection
private List<DetailBean> list;
DetailBean parameter description
// duration of a sedentary type
private int d;
//0 static (less than 30 steps in 1 minute), 1 trigger sedentary, 2 movement
(more than 30 steps)
private int t;
```

## Sync heart rate data

```
nowTime current time zone * 3600 + unix second value of current time
Sync yesterday: nowTime-(24*3600)*1,
Sync the day before yesterday: nowTime-(24*3600)*2
Data can be synchronized for up to three days
val time = (getTimeZone() * 3600).toInt()
val nowTime = date.unixTimestamp + time

CommandHandle.getInstance().executeReqCmd(
    ReadHeartRateReq(nowTime),
    ICommandResponse<ReadHeartRateRsp> {

})

ReadHeartRateRsp parameter description
//nothing yet
private int size = 0;
//nothing yet
private int index = 0;
// unix second value of heart rate data
private int mUtcTime;
//The heart rate data array is one point every 5 minutes, the data subscript *5
is equal to the number of minutes of the day
private byte[] mHeartRateArray;
private boolean endFlag = false;
```

## Synchronized blood pressure function

```
//Synchronized automatic blood pressure, measured once an hour
CommandHandle.getInstance()
    .executeReqCmd(SimpleKeyReq(Constants.CMD_BP_TIMING_MONITOR_DATA),
ICommandResponse<BpDataRsp> {})
BpDataEntity parameter description
//year
private int year;
```

```

//moon
private int month;
//day
private int day;
private int timeDelay;
private ArrayList<BpValue> bpValues;

BpValue parameter description
//The minute of the day, usually the whole hour
int timeMinute;
//measured heart rate value
int value;

Get the blood pressure value calculated from the measured value
//The heart rate value returned by the hr callback, age is the age of the user
val sbp= CalcBloodPressureByHeart.cal_sbp(hr, age) (systolic blood pressure)
//sbp heart rate calculated value
val dbp=CalcBloodPressureByHeart.cal_dbp(sbp) (diastolic pressure)

        //Confirm blood pressure synchronization, call after receiving the
        callback, the watch will delete the records that have been synchronized
        CommandHandle.getInstance().executeReqCmd(BpReadConformReq(true),null)

        //synchronize manual blood pressure
        CommandHandle.getInstance()
            .executeReqCmd(ReadPressureReq(0),
ICommandResponse<ReadBlePressureRsp> {}

        ReadBlePressureRsp.getValueList() parameter description
BlePressure parameter descriptionw
//time seconds value
public long time;
//(Diastolic pressure)
public int dbp;
//(systolic blood pressure)
public int sbp;

```

## Synchronized blood oxygen function

```

        LargeDataHandler.getInstance().syncBloodOxygenWithCallback(new
IBloodOxygenCallback() {
            @Override
            public void readBloodOxygen(List<BloodOxygenEntity> data) {

            }
});

```

BloodOxygenEntity parameter description

```

//data date
private String dateStr;
//Data minimum value array, one data per hour, a total of 24
private List<Integer> minArray;
//Data maximum value array, one data per hour, a total of 24
private List<Integer> maxArray;
//Data value at 0:00 on a certain day

```

```
private long unix_time;
```

## Synchronized Time Interval blood oxygen function

```
//One data point per minute, 1,440 data points per day
//dayIndex:Index of the day (0 = today, 1 = yesterday, and so on)
LargeDataHandler.getInstance().syncIntervalBloodOxygenWithCallback(dayIndex,new
IBloodOxygenCallback() {
    @Override
    public void readIntervalBloodOxygen(IntervalBloodOxygenEntity data) {

    }
});
```

IntervalBloodOxygenEntity parameter description

```
//dayIndex
private int dayIndex;
//interval
private int interval;
//Data value array, one data per min
private List<Integer> array;
```

## Synchronized Time Interval temperature function

```
//One data point per minute, 1,440 data points per day
//dayIndex:Index of the day (0 = today, 1 = yesterday, and so on)
LargeDataHandler.getInstance().syncIntervalTemperatureWithCallback(dayIndex,new
IIntervalTemperatureCallback() {
    @Override
    public void readIntervalTemperature(IntervalTemperatureEntity data) {

    }
});
```

IntervalTemperatureEntity parameter description

```
//dayIndex
private int dayIndex;
//interval
private int interval;
//Data value array, one data per min
private List<Float> array;
```

## Synchronized pressure function

```
//offset 0-6
//Synchronization pressure 7 days 0 only synchronizes today 1 yesterday 2
synchronizes the day before yesterday .... supports up to 7 days
CommandHandle.getInstance()
    .executeReqCmd(PressureReq(offset),
```



```

        ICommandResponse<PressureRsp> {

    })

    //PressureRsp Parameter Description
    //Pressure data, the return value divided by 10 is the value displayed by the
    APP, and one data is generated every half hour.
    private byte[] pressureArray;
    //Stress test interval
    private int range=30;
    //Pressure data date
    private DateUtil today;

```

## Synchronized hrv function

```

    //offset 0-6
    //Synchronization hrv 7 days 0 only synchronizes today 1 yesterday 2 synchronizes
    the day before yesterday .... supports up to 7 days
    CommandHandle.getInstance()
        .executeReqCmd(HRVReq(0),
            ICommandResponse<HRVRsp> {

    })

    //HRVRsp Parameter Description
    //hrv data, the return value divided by 10 is the value displayed by the APP, and
    one data is generated every half hour.
    private byte[] pressureArray;
    //hrv test interval
    private int range=30;
    //hrv data date
    private DateUtil today;

```

## Synchronized skin(body) temperature function

```

    fun registerTempCallback(){
        //Register a callback
        FileHandle.getInstance().clearCallback()
        FileHandle.getInstance().registerCallback(Callback())
        FileHandle.getInstance().initRegister()
    }

    //Synchronized automatic skin temperature
    fun syncAutoTemperature(){
        ////Synchronize automatic body temperature for 3 days 0 Synchronize only
        today 1 Synchronize today and yesterday 2 Synchronize today, yesterday and the
        day before yesterday .... Support up to 7 days
        FileHandle.getInstance().startObtainTemperatureSeries(2)
    }

    //Synchronous manual skin temperature
    fun syncManual(){
        //Synchronize automatic body temperature for 3 days 0 Synchronize only today
        1 Synchronize today and yesterday 2 Synchronize today, yesterday and the day
        before yesterday .... Support up to 7 days
        FileHandle.getInstance().startObtainTemperatureOnce(0)
    }

```

```

}

inner class callback : SimpleCallback() {
    //Continuous temperature callback
    override fun onUpdateTemperature(data: TemperatureEntity) {

    }
    //Single temperature callback
    override fun onUpdateTemperatureList(array:
MutableList<TemperatureOnceEntity>) {

    }
}
}

TemperatureEntity Parameter Description
public int mIndex; represents today 1 represents yesterday ...6
public int mTimeSpan; Test interval
public float[] mValues; Temperature array cursor*mTimeSpan = measurement time of
the day, in minutes
public float[] mSideValues; Side temperature array cursor*mTimeSpan =
measurement time in a day, in minutes

TemperatureOnceEntity Parameter Description
public long mTime; Measurement timestamp, seconds
public float mValue;Measurement value

```

## Synchronous training records

```

val syncSport= SportPlusHandle()
syncSport.timeFormat="yyyy-MM-dd HH:mm"
syncSport.syncSportPlus(object :
BaseTransientBottomBar.BaseCallback<MutableList<SportPlusEntity>>() {
    public override fun result(
        errorCode: Int,
        sportPlusEntities: MutableList<SportPlusEntity>

    ) {
    }
})

//Obtain all the data from the specified time until the present.0 is mean all.
//Note:The equipment only stores the latest ten pieces of data.
syncSport.cmdSummary(0)

SportPlusEntity
/**
 * Movement type index
 */
public int mSportType;
/**
 * Start time seconds
 */
public int mStartTime;
/**

```

```

    * Movement duration seconds
    */
    public int mDuration;
    /**
     * Movement mileage meters
     */
    public int mDistance;
    /**
     * Calories small calories
     */
    public float mCalories;
    /**
     * Average speed cm/s
     */
    public int mSpeedAvg;
    /**
     * Maximum speed cm/s
     */
    public int mSpeedMax;
    /**
     * Average heart rate beats/minute
     */
    public int mRateAvg;
    /**
     * Minimum heart rate beats/minute
     */
    public int mRateMin;
    /**
     * Maximum heart rate beats/minute
     */
    public int mRateMax;
    /**
     * Average altitude cm
     */
    public int mElevation;
    /**
     * Cumulative climbing cm
     */
    public int mUphill;
    /**
     * Cumulative downhill cm
     */
    public int mDownhill;
    /**
     * Average cadence steps/minute
     */
    public int mStepRate;
    /**
     * Number of exercises times
     */
    public int mSportCount;

    public int steps;

    public List<SportLocation> mLocations = new ArrayList<>();

```

```

SportLocation {

/**
 * Real-time heart rate beats/minute
 */
    public int mRateReal;
}

typedef NS_ENUM(NSInteger, OdmSportPlusExerciseModelType) {
    OdmSportPlusExerciseModelTypeGpsRun = 1, //gps running
    OdmSportPlusExerciseModelTypeGpsBike = 2, //gps cycling
    OdmSportPlusExerciseModelTypeGpsWalk = 3, //gps walking
    OdmSportPlusExerciseModelTypeIndoorRun = 4, //bracelet walking
    OdmSportPlusExerciseModelTypeRopeSkipping = 5, //bracelet rope
    skipping
    OdmSportPlusExerciseModelTypeSwimming = 6, //bracelet swimming
    OdmSportPlusExerciseModelTypeRun = 7, //bracelet running-outdoor
    OdmSportPlusExerciseModelTypeWalk = 8, //Bracelet hiking
    OdmSportPlusExerciseModelTypeBike = 9, //Bracelet cycling
    OdmSportPlusExerciseModelTypeExercise = 10, //Others
    OdmSportPlusExerciseModelTypeSwing = 11, //Bracelet swinging
    OdmSportPlusExerciseModelTypeClimb = 20, //Bracelet climbing
    OdmSportPlusExerciseModelTypeBadminton = 21, //Bracelet badminton
    OdmSportPlusExerciseModelTypeYoga = 22, //Bracelet yoga
    OdmSportPlusExerciseModelTypeAerobics = 23, //Bracelet aerobics
    OdmSportPlusExerciseModelTypeSpinningBike = 24, //Bracelet spinning
    bike
    OdmSportPlusExerciseModelTypeKayaking = 25, // Bracelet kayaking
    OdmSportPlusExerciseModelTypeEllipticalMachine = 26, // Bracelet
    elliptical machine
    OdmSportPlusExerciseModelTypeRowingMachine = 27, // Bracelet rowing
    machine
    OdmSportPlusExerciseModelTypePingpong = 28, // Bracelet table tennis
    OdmSportPlusExerciseModelTypeTennis = 29, // Bracelet tennis
    OdmSportPlusExerciseModelTypeGolf = 30, // Bracelet golf
    OdmSportPlusExerciseModelTypeBasketball = 31, // Bracelet basketball
    OdmSportPlusExerciseModelTypeFootball = 32, // Bracelet football
    OdmSportPlusExerciseModelTypeVolleyball = 33, // Bracelet volleyball
    OdmSportPlusExerciseModelTypeRockClimbing = 34, // Bracelet rock
    climbing
    OdmSportPlusExerciseModelTypeDance = 35, // Bracelet dance
    OdmSportPlusExerciseModelTypeRollerSkating = 36, // Bracelet roller
    skating

    OdmSportPlusExerciseModelTypeTreadmill = 40, // Running-treadmill
    OdmSportPlusExerciseModelTypeIndoorWalking = 41, // Running-indoor
    walking
    OdmSportPlusExerciseModelTypeTrailRunning = 42, // Running-cross-
    country running
    OdmSportPlusExerciseModelTypeRaceWalk = 43, // Running-race walking
    OdmSportPlusExerciseModelTypePlaygroundRunning = 44, // Running -
    playground running
    OdmSportPlusExerciseModelTypeFatLossRunning = 45, // Running - fat
    loss running

```

```

OdmSportPlusExerciseModelTypeOutdoorCycling = 50, // Cycling -
outdoor cycling
OdmSportPlusExerciseModelTypeIndoorCycling = 51, // Cycling - indoor
cycling
OdmSportPlusExerciseModelTypeMountainBiking = 52, // Cycling -
mountain biking
OdmSportPlusExerciseModelTypeBMX = 53, // Cycling - BMX

OdmSportPlusExerciseModelTypeSwimmingPool = 55, // Swimming -
swimming pool swimming
OdmSportPlusExerciseModelTypeOutdoorSwimming = 56, // Swimming -
outdoor swimming
OdmSportPlusExerciseModelTypeFinSwimming = 57, // Swimming-Fin
Swimming
OdmSportPlusExerciseModelTypeSynchronizedSwimming = 58, // Swimming-
Synchronized Swimming

OdmSportPlusExerciseModelTypeOutdoorHiking = 60, // Outdoor Sports-
Outdoor Hiking
OdmSportPlusExerciseModelTypeOrienteering = 61, // Outdoor Sports-
Orienteering
OdmSportPlusExerciseModelTypeFishing = 62, // Outdoor Sports-Fishing
OdmSportPlusExerciseModelTypeHunt= 63, // Outdoor Sports-Hunting
OdmSportPlusExerciseModelTypeSkateboard = 64, // Outdoor Sports-
skateboarding
OdmSportPlusExerciseModelTypeParkour = 65, // Outdoor Sports-Parkour
OdmSportPlusExerciseModelTypeATV = 66, // Outdoor sports - ATV
OdmSportPlusExerciseModelTypeMotocross = 67, // Outdoor sports -
Motocross
OdmSportPlusExerciseModelTypeRacing = 68, // Outdoor sports - Racing
OdmSportPlusExerciseModelTypeHandCrank = 69,
// Outdoor sports - hand bike
OdmSportPlusExerciseModelTypeMarathon = 70, // Outdoor sports -
marathon
OdmSportPlusExerciseModelTypeObstacleCourse = 71, // Outdoor sports -
obstacle course

OdmSportPlusExerciseModelTypeStairClimber = 80, // Indoor sports -
stair climbing machine
OdmSportPlusExerciseModelTypeStairStepper = 81, // Indoor sports -
stepper
OdmSportPlusExerciseModelTypeMixedAerobic = 82, // Indoor sports -
mixed aerobics
OdmSportPlusExerciseModelTypeKickboxing = 83, // Indoor sports -
kickboxing
OdmSportPlusExerciseModelTypeCoreTraining = 84, // Indoor sports -
core training
OdmSportPlusExerciseModelTypeCrossTraining = 85, // Indoor sports -
cross training
OdmSportPlusExerciseModelTypeIndoorFitness = 86, // Indoor sports -
indoor fitness
OdmSportPlusExerciseModelTypeGroupGymnastics = 87, // Indoor sports -
group gymnastics
OdmSportPlusExerciseModelTypeStrengthTraining = 88, // Indoor sports
- strength training

```

```

OdmSportPlusExerciseModelTypeGapTraining = 89, // Indoor sports - gap
training
OdmSportPlusExerciseModelTypeFreeTraining = 90, // Indoor sports -
free training
OdmSportPlusExerciseModelTypeFlexibilityTraining = 91, // Indoor
sports - flexibility training
OdmSportPlusExerciseModelTypeGymnastics = 92, // Indoor sports -
gymnastics
OdmSportPlusExerciseModelTypeStretch = 93, // Indoor exercise -
stretching
OdmSportPlusExerciseModelTypePilates = 94, // Indoor exercise -
Pilates
OdmSportPlusExerciseModelTypeHorizontalBar = 95, // Indoor exercise -
horizontal bar
OdmSportPlusExerciseModelTypeParallelBars = 96, // Indoor exercise -
parallel bars
OdmSportPlusExerciseModelTypeBattleRope = 97, // Indoor exercise -
battle rope
OdmSportPlusExerciseModelTypeFitness = 98, // Indoor exercise -
fitness
OdmSportPlusExerciseModelTypeBalanceTraining = 99, // Indoor exercise
- balance training
OdmSportPlusExerciseModelTypeStepTraining = 100, // Indoor exercise -
step training

OdmSportPlusExerciseModelTypeSquareDance = 110, // Dance sports-
square dance
OdmSportPlusExerciseModelTypeBallroomDancing = 111, // Dance sports-
social dance
OdmSportPlusExerciseModelTypeBellyDance = 112, // Dance sports-belly
dance
OdmSportPlusExerciseModelTypeBallet = 113, // Dance sports-ballet
OdmSportPlusExerciseModelTypeStreetDance = 114, // Dance sports-
street dance
OdmSportPlusExerciseModelTypeZumba = 115, // Dance sports-Zumba
OdmSportPlusExerciseModelTypeLatinDance = 116, // Dance sports-Latin
dance
OdmSportPlusExerciseModelTypeLatinJazz = 117, // Dance sports-Jazz
dance
OdmSportPlusExerciseModelTypeHipHopDance = 118, // Dance sports -
hip-hop dance
OdmSportPlusExerciseModelTypePoleDancing = 119, // Dance sports -
pole dance
OdmSportPlusExerciseModelTypeBreakDance = 120, // Dance sports -
break dance
OdmSportPlusExerciseModelTypeFolkDance = 121, // Dance sports - folk
dance
OdmSportPlusExerciseModelTypeNewDance = 122, // Dance sports - music
dance
OdmSportPlusExerciseModelTypeModernDance = 123, // Dance sports -
modern dance
OdmSportPlusExerciseModelTypeDisco = 124, // Dance sports - disco
OdmSportPlusExerciseModelTypeTapDance = 125, // Dance sports - tap
dance
OdmSportPlusExerciseModelTypeOtherDance = 126, // Dance sports -
other dances

```

```

OdmSportPlusExerciseModelTypeBoxing = 130, // Combat sports - boxing
OdmSportPlusExerciseModelTypeWrestling = 131, // Combat sports -
wrestling
OdmSportPlusExerciseModelTypeMartialArts= 132, // Combat sports -
martial arts
OdmSportPlusExerciseModelTypeTaiChi= 133, // Combat sports - Tai Chi
OdmSportPlusExerciseModelTypeMuayThai= 134, // Combat sports - Muay
Thai
OdmSportPlusExerciseModelTypeJudo= 135, // Combat sports - Judo
OdmSportPlusExerciseModelTypeTaekwondo= 136, // Combat Sports -
Taekwondo
OdmSportPlusExerciseModelTypeKarate= 137, // Combat Sports - Karate
OdmSportPlusExerciseModelTypeFreeSparring= 138, // Combat Sports -
Free Fighting
OdmSportPlusExerciseModelTypeSwordsmanship= 139, // Combat Sports -
Swordsmanship
OdmSportPlusExerciseModelTypeJujitsu= 140, // Combat Sports - Jiu-
Jitsu
OdmSportPlusExerciseModelTypeFencing= 141, // Combat Sports - Fencing
OdmSportPlusExerciseModelTypeKendo= 142, // Combat Sports - Kendo

OdmSportPlusExerciseModelTypeBeachFootball=150, //!<150Ball Sports -
Beach Soccer
OdmSportPlusExerciseModelTypeBeachVolleyball=151, //!<151Ball Sports-
Beach volleyball
OdmSportPlusExerciseModelTypeBaseball=152, //!<152Ball Sports-
Baseball
OdmSportPlusExerciseModelTypeSoftball=153, //!<153Ball Sports-
Softball
OdmSportPlusExerciseModelTypeNewFootball=154, //!<154Ball Sports-
Rugby
OdmSportPlusExerciseModelTypeHockey=155, //!<155Ball Sports-Hockey
OdmSportPlusExerciseModelTypeSquash=156, //!<156Ball Sports-Squash
OdmSportPlusExerciseModelTypeDoorKick=157, //!<157Ball Sports-
Gateball
OdmSportPlusExercis
eModelTypeCricket=158, //!<158Ball Sports-Cricket
OdmSportPlusExerciseModelTypeHandball=159, //!<159Ball Sports-
Handball
OdmSportPlusExerciseModelTypeBowling=160, //!<160Ball Sports-Bowling
OdmSportPlusExerciseModelTypePolo=161, //!<161Ball Sports-Polo
OdmSportPlusExerciseModelTypeRacquetball=162, //!<162Ball sports-
Racquetball
OdmSportPlusExerciseModelTypeBilliards=163, //!<163Ball Sports-Table
Tennis
OdmSportPlusExerciseModelTypeTakraw=164, //!<164Ball Sports-Sepak
Rattan Ball
OdmSportPlusExerciseModelTypeDodgeBall=165, //!<165Ball Sports-
Dodgeball
OdmSportPlusExerciseModelTypeWaterPolo=166, //!<166Ball Sports-Water
Polo
OdmSportPlusExerciseModelTypePuck=167, //!<167Ball Sports-Ice Hockey
OdmSportPlusExerciseModelTypeShuttlecock=168, //!<168Ball Sports-
Shuttlecock

```

OdmSportPlusExerciseModelTypeIndoorSoccer=169, //!<169Ball Sports-Indoor Football

OdmSportPlusExerciseModelTypeSandbag=170, //!<170Ball Sports-Sandbag Ball

OdmSportPlusExerciseModelTypeBocce=171, //!<171Ball Sports-Floor Bocce

OdmSportPlusExerciseModelTypeJaiBall=172, //!<172Ball Sports-Paiai

OdmSportPlusExerciseModelTypeFloorBall=173, //!<173Ball Sports-Floorball

OdmSportPlusExerciseModelTypeAustralianRulesFootball=174, //!<174Ball Sports-Australian Rules Football

OdmSportPlusExerciseModelTypePickering=175, //!<175Ball Sports-Pickering

OdmSportPlusExerciseModelTypeOutdoorBoating=180, //!<180Water Sports-Outdoor Rowing

OdmSportPlusExerciseModelTypeSailing=181, //!<181Water Sports-Sailing

OdmSportPlusExerciseModelTypeDragonBoat=182, //!<182 Water Sports-Dragon Boat

OdmSportPlusExerciseModelTypeSurf=183, //!<183 Water Sports-Surfing

OdmSportPlusExerciseModelTypeKitesurfing=184, //!<184 Water Sports-Kite Surfing

OdmSportPlusExerciseModelTypePaddling=185, //!<185 Water Sports-Paddling

OdmSportPlusExerciseModelTypePaddleboard=186, //!<186 water Sports-Paddleboard Surfing

OdmSportPlusExerciseModelTypeIndoorSurfing=187, //!<187 water Sports-Indoor Surfing

OdmSportPlusExerciseModelTypeDrifting=188, //!<188 water Sports-Rafting

OdmSportPlusExerciseModelTypeSnorkeling=189, //!<189 Water Sports - Snorkeling

OdmSportPlusExerciseModelTypeSkis=190, //!<190 Snow Sports - Double Skiing

OdmSportPlusExerciseModelTypeSnowboard=191, //!<191 Snow Sports - Singleboarding

OdmSportPlusExerciseModelTypeAlpineSkiing=192, //!<192 Snow Sports - Alpine Skiing

OdmSportPlusExerciseModelTypeCrossCountrySkiing=193, //!<193 Snow Sports - Cross-Country Skiing

OdmSportPlusExerciseModelTypeSkiOrientserIng=194, //!<194 Snow Sports - Orienteering

OdmSportPlusExerciseModelTypeBiathlon=195, //!<195 Snow Sports - Biathlon

OdmSportPlusExerciseModelTypeOutdoorskating=196, //!<196Snow and ice sports-outdoor skating

OdmSportPlusExerciseModelTypeIndoorskating=197, //!<197Snow and ice sports-indoor skating

OdmSportPlusExerciseModelTypeCurling=198, //!<198Snow and ice sports-curling

OdmSportPlusExerciseModelTypeBobsleigh=199, //!<199Snow and ice sports-bobsleigh

OdmSportPlusExerciseModelTypeSled=200, //!<200Snow and ice sports-sled



```

        OdmSportPlusExerciseModelTypeSnowmobile=201, //!<201Snow and ice
sports-snowmobile
        OdmSportPlusExerciseModelTypeSnowshoeing=202, //!<202 Ice and Snow
Sports-Snowshoe Hiking

        OdmSportPlusExerciseModelTypeHulaHoop=210, //!<210 Leisure Sports-
Hula Hoop
        OdmSportPlusExerciseModelTypeFrisbee=211, //!<211 Leisure Sports-
Frisbee
        OdmSportPlusExerciseModelTypeDarts=212, //!<212 Leisure Sports-Darts
        OdmSportPlusExerciseModelTypeFlyAKite=213, //!<213 Leisure Sports-
Kite Flying
        OdmSportPlusExerciseModelTypeTugOfWar=214, //!<214 Leisure Sports-Tug
of War
        OdmSportPlusExerciseModelTypeEsports=215, //!<215 Leisure Sports-
Esports
        OdmSportPlusExerciseModelTypeStroller=216, //!<216 Leisure Sports-
Walking Machine
        OdmSportPlusExerciseModelTypeNewSwing=217, //!<217 Leisure Sports-
Swing
        OdmSportPlusExerciseModelTypeShuffleboard=218, //!<218 Leisure
Sports-Shuffleboard
        OdmSportPlusExerciseModelTypeTableSoccer=219, //!<219 Leisure Sports-
Table Football
        OdmSportPlusExerciseModelTypeSomatosensoryGame=220, //!<220 Leisure
Sports-Somatics
        OdmSportPlusExerciseModelTypeBungeeJumping=221, //!<221 Leisure
Sports-Bungee Jumping
        OdmSportPlusExerciseModelTypeParachute=222, //!<222 Leisure Sports-
Skydiving
        OdmSportPlusExerciseModelTypeAnusara=223, //!<223 Leisure Sports-
Anusara
        OdmSportPlusExerciseModelTypeYinYoga=224, //!<224 Leisure Sports-Yin
Yoga
        OdmSportPlusExerciseModelTypePregnancyY
..

```

### 2.3.6 OTA upgrade function:

```

```java
    //dfu upgrade instance
    val fuHandle= DfuHandle.getInstance()
    //initialize callback
    dfuHandle.initCallback()
    //DFU file verification, path firmware file path
    if (dfuHandle.checkFile(path)) {
        dfuHandle.start(dfuOpResult)
    }
    //dfuOpResult callback description
    private val dfuOpResult: DfuHandle.IOpResult = object : DfuHandle.IOpResult {
        override fun onActionResult(type: Int, errCode: Int) {
            if (errCode == DfuHandle.RSP_OK) {
                when (type) {
                    1 -> dfuHandle.init()
                    2 -> dfuHandle.sendPacket()
                    3 -> dfuHandle.check()

```

```

        4 -> {
            //The upgrade is successful, wait for the device to
restart
            dfuHandle.endAndRelease()
        }
    }
} else {
    //Upgrade exception or failure
}
}

override fun onProgress(percent: Int) {
    // file upgrade progress
}
}

```

## 2.3.7 Manual measurement

```

//StartHeartRateRsp parameter description
private byte type; type 1: heart rate 2: blood pressure 3: blood oxygen
private byte errCode; measurement error code 0: normal 1: measurement failed 2:
measurement failed
private byte value; measurement value: heart rate or blood oxygen
private byte sbp; blood pressure sbp
private byte dbp; blood pressure dbph

Boolean stop true stop false start

// manual heart rate
BleOperateManager.getInstance().manualModeHeart(new
ICommandResponse<StartHeartRateRsp>() {
    @Override
    public void onDataResponse(StartHeartRateRsp resultEntity) {

    }
}, Boolean stop);
    //manual blood pressure
    BleOperateManager.getInstance().manualModeBP(new
ICommandResponse<StartHeartRateRsp>() {
    @Override
    public void onDataResponse(StartHeartRateRsp resultEntity) {

    }
}, Boolean stop);
    //manual blood oxygen
    BleOperateManager.getInstance().manualModeSpO2(new
ICommandResponse<StartHeartRateRsp>() {
    @Override
    public void onDataResponse(StartHeartRateRsp resultEntity) {

    }
}, Boolean stop);
    //manual pressure

```

```

        BleOperateManager.getInstance().manualModePressure(new
ICommandResponse<StartHeartRateRsp>() {
    @Override
    public void onDataResponse(StartHeartRateRsp resultEntity) {

    }
},Boolean stop);
        //manual hrv
        BleOperateManager.getInstance().manualModeHrv(new
ICommandResponse<StartHeartRateRsp>() {
    @Override
    public void onDataResponse(StartHeartRateRsp resultEntity) {

    }
},Boolean stop);
        //manual temperature
        BleOperateManager.getInstance().manualTemperature(new
ICommandResponse<StartHeartRateRsp>() {
    @Override
    public void onDataResponse(StartHeartRateRsp resultEntity) {
        //The temperature value obtained should be divided by 10 to get the
normal temperature.
    }
},Boolean stop);
        //One-click measurement

        BleOperateManager.getInstance().oneClickMeasurement(new
ICommandResponse<StopHeartRateRsp>() {
    @Override
    public void onDataResponse(StopHeartRateRsp resultEntity) {

    }
},Boolean stop);

        BleOperateManager.getInstance().manualModeHeartRateRawData(ICommandResponse<Stop
HeartRateRsp> response,int seconds,boolean stop) (new
ICommandResponse<StopHeartRateRsp>() {
    @Override
    public void onDataResponse(StopHeartRateRsp resultEntity) {

    }
},Boolean stop);

        BleOperateManager.getInstance().manualModeBloodOxygenRawData(ICommandResponse<St
opHeartRateRsp> response,int seconds,boolean stop) (new
ICommandResponse<StopHeartRateRsp>() {
    @Override
    public void onDataResponse(StopHeartRateRsp resultEntity) {

    }
},Boolean stop);

//StopHeartRateRsp params
private int sbp;

```

```

private int dbp;
private int rri;
private int heart;
private int hrv;
private int stress;
private int temperature;

private int ppgCount;
private int greenLightPpgL;
private int greenLightPpgH;
private int XL;
private int XH;
private int YL;
private int YH;
private int ZL;
private int ZH;
private int heartRate;
private int redLightPpgL;
private int redLightPpgH;
private int infraredPpgL;
private int infraredPpgH;
private int bloodOxygen;

```

## 2.3.8 Touch and gestures

```

    read
touch: true  touch  false:gestures

    CommandHandle.getInstance().executeReqCmd(TouchControlReq.getReadInstance(false)
,
        ICommandResponse<TouchControlResp> {
            appType:  0:close  1:music  2:video  3:muslim  4:ebook  5:camera  6:phone call
7:game  8:heart
            strength: Dynamics
        }

write:
appType:  0:close  1:music  2:video  3:muslim  4:ebook  5:camera  6:phone call  7:game
8:heart  9:touch event
touch: true  touch  false:gestures
Strength:1-10

    CommandHandle.getInstance().executeReqCmdNoCallback(TouchControlReq.getWriteInsta
nce(appType,false,currStrength))

    //Set RT11 Touch Function

    CommandHandle.getInstance().executeReqCmd(TouchControlReq.getWriteTpSleepInstanc
e(appType,1),null)

```

## 2.3.9 Changes in Ring measurement data are proactively reported to the APP

```
//Add listener
```

```
BleOperateManager.getInstance().addOutDeviceListener(ListenerKey.Heart,myDeviceNotifyListener)
```

ListenerKey Parameter Description

```
public class ListenerKey {  
    public static int Heart=1;           Heart  
    public static int BloodPressure=2;   Blood Pressure  
    public static int BloodOxygen=3;     Blood Oxygen  
    public static int Temperature=5;     Temperature  
    public static int SportRecord=7;     Sport Record  
    public static int All=7;             All  
}
```

Monitoring instructions

```
inner class MyDeviceNotifyListener : DeviceNotifyListener() {  
    override fun onDataResponse(resultEntity: DeviceNotifyRsp?) {  
        if (resultEntity!!.status == BaseRspCmd.RESULT_OK) {  
            BleOperateManager.getInstance().removeOthersListener()  
            when (resultEntity.dataType) {  
                1 -> {  
                    //Watch heart rate test changes  
                }  
                2 -> {  
                    //Watch blood pressure test changes  
                }  
                3 -> {  
                    //Watch blood oxygen test changes  
                }  
                4 -> {  
                    //Changes in watch step counting details  
                }  
                5 -> {  
                    //Watch body temperature changes on the day  
                }  
                7 -> {  
                    //Generate new exercise records  
                }  
                0x0c -> {  
                    //730c5701000000000000000000000000d7  
                    val charging = BLEDataFormatUtils.bytes2Int(  
                        byteArrayOf(  
                            resultEntity.loadData[2]  
                        )  
                    )  
                    if (charging > 0) {  
                        //charging  
                    } else {  
                        val battery = BLEDataFormatUtils.bytes2Int(  
                            byteArrayOf(  
                                resultEntity.loadData[2]  
                            )  
                        )  
                    }  
                }  
            }  
        }  
    }  
}
```

```

                                resultEntity.loadData[1]
                            )
                        )
                        //battery power
                    }
                }

0x2d ->{
    //The custom function button is triggered
    val event = BLEDataFormatUtils.bytes2Int(
        byteArrayOf(
            resultEntity.loadData[1]
        )
    )
    event:0 null 1:decline 2: Slide up 3: Single 4: Long press
}

0x12 -> {
    //7312 00005200025100003c0000000066
    AwLog.i(Author.HeZhiYuan,
ByteUtil.bytesToString(resultEntity.loadData))

    val step = BLEDataFormatUtils.bytes2Int(
        byteArrayOf(
            resultEntity.loadData[1],
            resultEntity.loadData[2],
            resultEntity.loadData[3]
        )
    )

    val calorie = BLEDataFormatUtils.bytes2Int(
        byteArrayOf(
            resultEntity.loadData[4],
            resultEntity.loadData[5],
            resultEntity.loadData[6]
        )
    )

    val distance = BLEDataFormatUtils.bytes2Int(
        byteArrayOf(
            resultEntity.loadData[7],
            resultEntity.loadData[8],
            resultEntity.loadData[9]
        )
    )

    deviceNotification(step, distance, calorie)
}

0x3A -> { //must add Heart rate tooLow/High Reminder
Listener:ListenerKey.Heart

    //BleOperateManager.getInstance()
    //.addOutDeviceListener(ListenerKey.Heart,
myDeviceNotifyListener)

    //Heart Rate too low/high Reminder
    val type = BLEDataFormatUtils.bytes2Int(

```

```

        byteArrayOf(
            resultEntity.loadData[1]
        )
    )
    if (type == 1) { //too low

    } else if (type == 2) { //too high

    }
    //Heart Rate value
    val value = BLEDataFormatUtils.bytes2Int(
        byteArrayOf(
            resultEntity.loadData[2]
        )
    )
    }
    }
    }
    }
}

//Remove the listener. It must be removed after registration, otherwise multiple
callbacks will appear.
BLEOperateManager.getInstance().removeNotifyListener(ListenerKey.Heart)
//Remove all listeners
BLEOperateManager.getInstance().removeNotifyListener(ListenerKey.All)

```

## 2.3.10 APP opens exercise type

```

// status 1 Start movement 2 Pause 3 Continue 4 End 6 Movement start timestamp
//Sport type 4 walking 5 Jumping rope 7 Running 8 Hiking 9 Cycling 10 Other
sports 20 Hiking 21 Badminton
22 Yoga 23 Aerobics 24 Spinning 25 Kayaking 26 Elliptical machine 27 Rowing
machine 28 Table tennis 29 Tennis
30 Golf 31 Basketball 32 Football 33 volleyball 34 Rock climbing 35 Dance 36
Roller skating 60 Outdoor hiking
CommandHandle.getInstance().executeReqCmd(
    PhoneSportReq.getSportStatus(
        1, sportType
    ), gpsResponse
)

private var gpsResponse: ICommandResponse<AppSportRsp> =
    ICommandResponse<AppSportRsp> { resultEntity ->
        AwLog.i(Author.HeZhiYuan, resultEntity)
        if (resultEntity != null) {
            when (resultEntity.gpsStatus) {
                6 -> {
                    //Exercise start time (Unit second)
                }

                2 -> {
                    //Exercise pause
                }
            }
        }
    }

```

```

3 -> {
    // //Exercise continues
}

4 -> {
    //Exercise end
}
0x25 -> {
    //Muslim ring click real-time data
    //732500000013000700000000000000b2
    val count = BLEDataFormatUtils.bytes2Int(
        byteArrayOf(
            resultEntity.loadData[1],
            resultEntity.loadData[2],
            resultEntity.loadData[3],
            resultEntity.loadData[4],
        )
    )
}
}

//Report data during exercise
//Add motion data reporting and monitoring
BLEoperateManager.getInstance().addSportDeviceListener(0x78,
myDeviceSportNotifyListener)
//Remove sports data reporting monitoring
BLEoperateManager.getInstance().removeSportDeviceListener(0x78)
//Listening to inner classes
inner class MyDeviceNotifyListener : DeviceSportNotifyListener() {
    override fun onDataResponse(resultEntity: DeviceNotifyRsp?) {
        super.onDataResponse(resultEntity)

        if (resultEntity!!.status == BaseRspCmd.RESULT_OK) {
            //Movement duration, unit seconds
            val sportTime = bytes2Int(
                byteArrayOf(
                    resultEntity.loadData[2],
                    resultEntity.loadData[3]
                )
            )
            //Exercise real-time heart rate
            val heart = bytes2Int(
                byteArrayOf(
                    resultEntity.loadData[4]
                )
            )

            //The number of steps generated during exercise will only have
            data when the exercise type is 4, 7, or 8, otherwise it will be 0
            val step = bytes2Int(
                byteArrayOf(
                    resultEntity.loadData[5],
                    resultEntity.loadData[6],
                    resultEntity.loadData[7]
                )
            )
        }
    }
}

```



```

        )
    )
    //The distance generated during exercise will only have data when
    the exercise type is 4, 7, or 8, and the others will be 0.
    val distance = bytes2Int(
        byteArrayOf(
            resultEntity.loadData[8],
            resultEntity.loadData[9],
            resultEntity.loadData[10]
        )
    )
    //Calories generated during exercise
    val calorie = bytes2Int(
        byteArrayOf(
            resultEntity.loadData[11],
            resultEntity.loadData[12],
            resultEntity.loadData[13]
        )
    )

    //Error status during exercise
    val status = bytes2Int(
        byteArrayOf(
            resultEntity.loadData[1]
        )
    )
    val sportType = BLEDataFormatUtils.bytes2Int(
        byteArrayOf(
            resultEntity.loadData[0]
        )
    )
    if (status == 0x03) {
        //Not wearing
    }
}

}

/**
 * Convert the byte array to int type, with the high byte of the array first
 *
 * @param data
 * @return
 */
public static int bytes2Int(byte[] data) {
    int length = data.length;
    int res = 0;
    for (int i = 0; i < length; i++) {
        res |= (data[i] & 0xFF) << (8 * (length - 1 - i));
    }
    return res;
}

```

## Muslim Data Synchronization

//dayOffset 0: Today 1: Yesterday 2: The day before yesterday, supports synchronization for up to 7 days

```
CommandHandle.getInstance()
    .executeReqCmd(
        MuslimReq(0),
        ICommandResponse<MuslimRsp> {
    })
```

```
//MuslimRsp
private int size = 0;
private int index = 0;
private byte[] pressureArray; //Muslim data one data point per hour
private boolean endFlag = false;
private int range=60;
private DateUtil today;
private int offset=-1;
```

Real-time data reporting, ring single reporting, please refer to 2.3.9 Type 0x25

###

## Setting Ring user Id

```
//read
CommandHandle.getInstance()
    .executeReqCmd(
        PhoneIdReq.getReadInstance(),
object : ICommandResponse<PhoneIdRsp> {
    override fun onDataResponse(phoneIdRsp: PhoneIdRsp) {
        //userId:phoneIdRsp.getUserId()
    }
})
//write
CommandHandle.getInstance()
    .executeReqCmd(
        PhoneIdReq.getWriteInstance("123456789123"),
object : ICommandResponse<PhoneIdRsp> {
    override fun onDataResponse(phoneIdRsp: PhoneIdRsp) {

    }
})
```

## 设置个人信息

```
CommandHandle.getInstance().executeReqCmd(  
    TimeFormatReq.getWriteInstance(  
        12,24小时切换, 0为24 1为12,  
        公英, 0公 1英,  
        男 true,女 false, age, height, weight, 120, 90, 160  
    ), ICommandResponse<TimeFormatRsp>() {})
```